


User's Guide to ^1H , ^{13}C NMR

- Note: the default mouse button is the left button. Always use the left one unless told otherwise.
 - For help, see Dr. Jasperse, Hagen 411-I.
1. Login
 - a. double click on **jasperse** icon (or type in "jasperse")
(Research Users: You should use your boss's login.)
 - b. Password "**chem355**" (need underline in between)
 - c. double click on "**xwinnmr**"
 2. Sample Insertion/Lock/Tune
 - a. remove cap from spectrometer if needed, and then click **LIFT ON/OFF** key on upper lefthand corner of SCM keyboard (to right of computer) to lift lock sample
 - b. place your sample in sample holder, adjust position using depth gage, and place in spectrometer [DO NOT PUT SAMPLE IN WITHOUT THE SAMPLE HOLDER! YOU WILL BREAK YOUR SAMPLE AND WRECK THE INSTRUMENT!]
 - c. click **LIFT ON/OFF** key on SCM keyboard to lower sample
 3. Acquiring the spectrum
 - a. type "**xmac**"
 - a listing of suggestions will come up
 - at present, all of these assume CDCl_3 as solvent
 - instructors/researchers, these can be easily customized for you needs. See Jasperse
 - b. select the experiment of interest, normally **ah1-tune** or **c13-tune** or **c13-notune**
 - Note: if you are going to run both ^1H and ^{13}C on the same sample, you don't need to tune twice. Run "ah1-tune" first for hydrogen, then "c13-notune"
 - c. when asked for file name information, type your name into the name box and
 - d. type "**chem355**" into the "user" box.
 - e. click **SAVE**
 - f. click **COPY ALL** when the box comes up
 - The computer will now do everything for you: read in the correct parameters, lock, tune if specified, adjust the receiver gain, acquire the spectrum, phase the spectrum, and store the phasing information for automatic integration.
 - g. When "ns = 8 (or 128)" box pops up, hit return to accept default, or else enter something different for the number of scans
 - h. wait patiently until either an "**xmac:finished**" or "**abs finished**" message appears
 - Hopefully this whole process will take less than 5 minutes for proton or less than 8 minutes for a carbon spectrum. If 8 minutes have passed and still incomplete, see Jasperse...
 - i. click on the $\overline{\Delta}$ icon (upper right corner of the icon group) to adjust the vertical scale of the viewed spectrum. (For example, if your baseline looks flat, this will fix it!)
 4. Plot
 - a. click **PLOT** icon
 - b. hit return in response to any boxes that appear
 - c. To do **horizontal expansions**, **manual integrations** or **vertical expansions**, see instructions on page 2.
 5. Exiting:
 - a. Replace your sample with the default sample, as described in part 3.
 - b. Type "lock cdcl3"
 - c. type "**exit**"
 - d. Say OK if it asks you anything about closing things
 - e. put cursor outside of any boxes into the blue area, then press the *right* mouse button, click **Logout** and click **Yes**.

1. Plotting **Horizontal Expansions**


- Make sure that the cursor is somewhere on the spectrum.
- Click the **left mouse button**. You will now get a doubled arrow.
- Move the doubled arrow to the left end of the area you want to expand and click the **center mouse button** to define the **left boundary**.
- Move the doubled arrow to the right end of the area you want to expand and again click the **center mouse button** to define the **right boundary**.
- click **PLOT** icon
- To get back to the full expansion, click the  icon
- To get out of the "doubled arrow" mode, click the **left mouse button**

2. **Manual Integration**

- click **INTEGRATE**
- define the regions of interest (see horizontal expansion instructions above)
- click **RETURN** and save your integral regions
 - Sometimes you may wish to improve the "flatness" of the integral, or you may wish to assign calibration values of your own choosing. Do the following:
- put the arrow within the region of your integral, and click the **left mouse button**. The integral under consideration will then get a star by it.
- click **CALIBRATE** and respond accordingly
- adjust the **BIAS** in order to get the left side of the integral level
- then adjust the **SLOPE** to get the right side of the integral level

3. Reducing the noise in noisy, dilute ^{13}C spectra. "**Power Spectrum**" .

[Do not use for ^1H spectra!]

- After getting the normal spectrum, type "**ps**"
 - click on the  icon (upper right corner of the icon group) to adjust the vertical scale of the viewed spectrum.
 - click **PLOT** icon
- Note: The "ps" command can make plots look prettier, by de-emphasizing noise. It does so by squaring all signals, however, so it will also de-emphasize small peaks that are real. In addition, by changing the relative sizes of peaks, it is incompatible with integration.

4. **Vertical Expansions**

- Type "**cy**" and increase or decrease the default value as you see fit. Doubling will double the printed heights, tripling will triple the printed heights, etc.
 - At default, cy=14, and is set so that the tallest peak in the spectrum will be 14cm tall. Thus, if you are wanting to expand a peak that is too tall, you need to multiply the cy as needed.

5. **Manual Phasing**

- click **PHASE**
- Click **BIGGEST**
- Click **PHO**, and keep finger held down
- drag, to adjust phase of biggest, marked peak
- Click **PH1**, and keep finger held down
- drag, to adjust phase of peaks distant from biggest
- click **RETURN**
- type "abs" if you want integrations to be automatically printed as a result

6. **Printing Titles**

- | | |
|--|--|
| a. type " setti " (for "set title") | b. delete existing title and type in new one |
| c. click save | d. click quit |
| | e. type " title " |
| | f. choose yes |